



Software Supporting X-Cite 120PC / 120PC XL

Manufacturer	Software	Release	Driver Available
Carl Zeiss	AIM	4.2	yes
Carl Zeiss	Zen	-	yes
Improvision	Openlab	-	no - controlled by RS232 commands via "Automations"
Intelligent Imaging Innovations (3i)	Slidebook	post-March 2006	yes
Molecular Devices	Meta Imaging Series	4.5.2 and up	yes
Media Cybernetics	Scope-Pro (with Image-Pro 5.x)	5.x	yes
Nikon	NIS - Elements	-	yes
Olympus	SIS	-	yes
Scanalytics/BD BioSciences	IPLab	-	no - controlled by RS232 commands in user-defined script

Notes:

1. Drivers can be obtained from the software supplier.
2. RS-232 commands can be found in the X-Cite 120PC / 120PC XL user guide



Melb : (03) 9480 4999
Syd : (02) 9705 8059
Email: sales@scitech.com.au
www.scitech.com.au