



PRO 2300B Binocular or PRO 2300T Trinocular Microscopes

The Pro 2300 series are a high quality microscope with variable 6V 20 W Halogen illumination.
Coaxial Mechanical Stage.

5 position nosepiece

Coaxial fine and course focussing adjustment with friction and max height setting.

Abbe 1.25 condenser with iris and swing out low power lens, on rack and pinion and centerable.

Eyepieces WF 10x (FN 18)

Objectives

Plan achromat 4x 0.10

Plan achromat 10x 0.25

Plan achromat 40x 0.65

Plan Achromat HI 100x 1.25 oil immersion

Accessories:

Phase contrast Kits

Dark field condensers, or non immersion or oil immersion type

TV C mounts or Photo mounts