

SMZ 6745 series



SMZ 6745[B/T]-1



SMZ 6745[B/T]-3



SMZ 6745[B/T]-5

SMZ6745B-1 ZOOM

Stereomicroscope,
Zoom range 0.67 – 4.5x
Working Distance 100mm
Eyepiece WF10X./22mm
Binocular Tube Inclined 45 deg.

SMZ6745B-3

with Incident and Transmitted Illumination

SMZ6745T-3

Zoom stereoscopic microscope with Trinocular (Photo) Tube
Zoom range 0.67 – 4.5x
Eyepiece WF10x/22mm
with Incident and Transmitted Light Illumination.

SMZ6745T-5

Zoom stereomicroscope with Trinocular (Photo) Tube
Zoom range 0.67 – 4.5x
Eyepiece WF10x/22mm
with Universal Stand

Accessories:

SZ-CTV	C Mount
SZ-PH	Photo Tube
SZ-L1	6v 15w Light with Transformer
SZ-L2	Fluorescent 8 watt circular light
SZ-DK	Dark Field Condenser
SZ-GC	Gem Clamp
SZ-0.5X	Supplementary Lens
SZ-1.5X	Supplementary Lens
SZ-2X	Supplementary Lens