



See things you've never seen before.

## Dual Mode Fluorescence Module Specifications

**CytoViva™**

### Dimensions

Length = 129mm  
Width = 88mm  
Height = 144mm  
Base = 140mm x 140mm

### Weight

1.75kg ( 62oz.)

### Filter Compatibility

Standard 25mm excitation filters. Emission filter dimensions dependent on microscope requirements.

### Light Guide

Supports 5mm OD (3mm bundle) liquid light guide.

### Light Source

The CytoViva DMF module has been designed to be compatible with the Exfo X-Cite 120 fluorescent light source.

### Warranty

One (1) year from date of purchase.



**Dual Mode Fluorescence Module**  
pictured with CytoViva High Resolution Adapter



**SciTech**  
the imaging specialists

Melb : (03) 9480 4999  
Syd : (02) 9705 8059  
Email: [sales@scitech.com.au](mailto:sales@scitech.com.au)  
[www.scitech.com.au](http://www.scitech.com.au)